



ANTIBACTERIAL DRESSING

**MEDIHONEY<sup>®</sup> HCS**

## MEDIHONEY<sup>®</sup> in Hydrogel Dressing with Super Absorbent Polymer




### FEATURES

- 63% Antibacterial *Leptospermum* (Manuka) Honey
- Helps maintain moist environment conducive to wound healing
- For dry to moderately exuding wounds
- Available in both Non-Adhesive and Adhesive
- Adhesive dressing does not require secondary dressing

### INDICATIONS

- Diabetic foot ulcers
- Leg ulcers
- Pressure ulcers (I-IV)
- 1st and 2nd degree burns
- Donor site
- Traumatic wounds
- Surgical wounds

### MEDIHONEY<sup>®</sup> HCS dressings ability to handle fluid vs. other absorbent dressings



Fluid absorption to EN 13726-2:2002 standard, in relative humidity of 55%

[www.dermasciences.com](http://www.dermasciences.com)

Derma Sciences Europe Ltd., Suite 2, 1A Park Street, Maidenhead SL6 1SL, U.K., Tel: +44 (0) 1628 625 916 • Fax: +44 (0) 1628 671 112  
 Medihoney<sup>®</sup> is a registered trademark of Medihoney Pty Ltd. and licensed by Comvita New Zealand Ltd.



## MEDIHONEY® in Hydrogel Dressing with Super Absorbent Polymer



CODE	DESCRIPTION	PRODUCT DIMENSIONS	QTY/BOX	BOX/CASE	QTY/CASE
780	Non-Adhesive	6cm x 6cm	10	5	50
781	Non-Adhesive	11cm x 11cm	10	5	50
782	Adhesive	7.2cm x 7.2cm, outer 11cm x 11cm	10	5	50
783	Adhesive	11.5 cm x 11.5cm, outer 15cm x 15cm	10	5	50

### DIRECTIONS FOR USE

1. Protect vulnerable wound edges with a suitable barrier cream such as MEDIHONEY® Barrier Cream.
2. Use an appropriately sized MEDIHONEY® HCS dressing.
3. Open pouch and remove dressing. Apply MEDIHONEY® HCS over the wound bed to achieve the intimate contact with the entire wound bed. (You may overlap wound margins and intact skin). MEDIHONEY® HCS may also be cut or fenestrated around a tube site.
4. **Non-Adhesive MEDIHONEY® HCS** may be cut to a desired size. Use a thin film or island/composite to cover/secure the Non-Adhesive HCS to the wound bed.  
*Note: Adhesive MEDIHONEY® HCS does not need a secondary dressing.*
5. The wound may be cleaned between dressing changes in accordance with standard wound cleaning protocol.

**Frequency of change:** The **MEDIHONEY® HCS Adhesive and Non-Adhesive** Dressings require changing when the dressing has been diluted by wound exudate. The dressing may be left in place for up to 7 days depending on wound exudate. If maceration of the surrounding skin occurs more frequent changing of the dressing is required.

**Note: As noted above, there may be increased levels of exudate upon initial use of the dressing. In the case of increased levels of exudate, protect the peri-wound skin with MEDIHONEY® Barrier Cream.**



[www.dermasciences.com](http://www.dermasciences.com)

Derma Sciences Europe Ltd., Suite 2, 1A Park Street, Maidenhead SL6 1SL, U.K., Tel: +44 (0) 1628 625 916 • Fax: +44 (0) 1628 671 112  
Medihoney® is a registered trademark of Medihoney Pty Ltd. and licensed by Comvita New Zealand Ltd.