

ADHESIVE BARRIER WOUND DRESSING

COVADERM[®] PLUS

COVADERM[®] PLUS



- Ideal composite dressing for wounds of all types
- Replace gauze, tape and scissors with an easier solution
- Rounded corners make a difference
- Variety of dressing sizes to fit any size wound
- Easy to use

 DeRoyal[®]



ADHESIVE BARRIER WOUND DRESSING

COVADERM® PLUS

FEATURES & BENEFITS

- Multi-layered composite wound dressing
- Maintains moist environment with protection from outside containments
- Four layer construction:
 - Protective non-adherent contact layer
 - Soft pad for drainage absorption
 - Breathable, semi-occlusive polyurethane barrier film
 - Conformable adhesive-coated fabric tape
- Prevents drainage strike-through
- Barrier to external contamination
- Minimizes wound trauma
- Provides excellent moisture vapor transmission
- Ideal for many wound types
- Available with radial slit for tube sites
- Primary or secondary dressing

INDICATIONS

- Pressure sores
- Leg Ulcers
- Chronic secreting lesions
- I.V. sites
- Laparotomy sites
- Arthroscopic incisions
- Tube sites


Product #	Tape	Pad	Qty/Bx	Qty/Cs
46-411	2½" x 2½"	1½" x 1½"	10	50
46-401	4" x 4"	2½" x 2½"	10	50
46-402	6" x 6"	4" x 4"	10	50
46-403	6" x 8"	4" x 6"	10	50
46-415	4" x 14"	2" x 11"	n/a	100
46-408	4" Diameter	2½" Diameter	25	200

TO PLACE AN ORDER OR FOR MORE INFORMATION
 PLEASE CALL: **800-251-9864**
 865.938.7828 (OUTSIDE USA)
 800.327.0304 (FAX)

APPLICATION


CHRONIC WOUNDS:

Covaderm Plus is ideally sized for most chronic wounds from pressure sores and leg ulcers to other secreting lesions.




I.V. SITES:

Use Covaderm Plus to cover a discontinued I.V. Site.



SURGICAL WOUNDS:

Use Covaderm Plus for dressing laparoscopy sites after surgery. Use the larger Covaderm Plus sizes for incisions that cover a broad surface area. Such as incisions associated with c-sections, saphenous vein harvest, open heart surgery, etc.



ORTHOPEDICS:

The 2" x 2" Covaderm Plus is perfect for covering incisions made during arthroscopic knee surgery. Use the 4" x 4" to cover incisions made during ankle surgery and to provide a secure, comfortable covering over joint pins under orthopedic casts.



TUBE SITES:

Covaderm Plus is also available in a circular configuration with a radial slit for use around tube sites.



DeRoyal®

200 DeBusk Lane
 Powell, Tennessee 37849
www.deroyal.com