

The dual – purpose
Repose® pump ensures
 optimal performance



Patented “smart valve”

Patented “smart valve” technology automatically senses the optimum internal pressure, exhausts excess air, and prevents over inflation.

Audible click

An audible click confirms the correct pressure has been set and provides reassurance that Repose® is ready to use.

Dual purpose Packaging and Pump

This dual purpose device incorporates clear instruction for inflation, cleaning, deflation, repacking, and includes an owner identification label.

Integrated instructions for use

Detailed instructions for use including limitations, precautions and cleaning instructions are all conveniently located on the inside of the pump.

In - built safety

The locking button ensures that the pump does not open during storage or transportation and the Safety Grip built into the base provides a firm grip during pumping.

Transportation and Storage

Compact and lightweight tube for low cost storage and transportation.



Repose® is CLINICALLY EFFECTIVE. Many experts and users confirm this statement.

An overview of the most relevant clinical papers and technical tests:

- **Challenging the pressure sore paradigm.**
Price P et al. Journal of Wound Care. April Vol 8, No4, 1999.
- **Two clinical evaluations of the Repose® system**
Ann MacFarlane, Sue Sayer. Wounds UK 2006, Vol 2, N°3.
- **Clinical evaluation of the effectiveness of the Repose® System.**
Osterbrink J et al. March 2005 (Florida International University, Anesthesiology Nursing Program).
- **Het belang van zithouding en drukreducerende kussens in het ontstaan van drukletsels.**
Prof. Dr. Tom Defloor, Prof. Dr. Mieke Grypdonck. Dep. of Nursing Science, Univerisity of Gent, Belgium.
- **Antidecubituskussens, drukvermindering of toch niet?**
Prof. Dr. Tom Defloor, Prof. Dr. Mieke Grypdonck. Dep. of Nursing Science, Univerisity of Gent, Belgium.
- **Recognizing the feet as being at risk from pressure damage.**
Bale S et al. British Journal of Nursing, 2001, Vol 10, N° 20
- **Pressure Relief For Heels: An Effective Innovation.**
Melhuish JP et al. Wound Healing Research Unit. University of Wales College of Medicine, Heathe Park, Cardiff, UK.
- **Pressure Ulcer Prevalence Audit: What are the benefits of doing it?**
Anne Ballard Wilson. Tissue Viability Nurse Fife Acute Operational Division. Scotland
Poster bekroond met de 2de prijs tijdens de EPUAP in Berlin 31/08 - 02/09 2006
- **Repose: the cost-effective solution for prompt discharge of patients**
Hampton S. British Journal of Nursing, 2000, Vol 9, N°21.
- **The use of a new overlay mattress in patients with chronic pain: impact on sleep and self-reported pain.**
Price P et al. Clinical Rehabilitation 2003; 17:488-492.
- **Bestemming van de drukreducerende eigenschappen van Repose®.**
Berlin Cert Prüf- und Zertifizierungsstelle für Medizinprodukte GmbH, 2006-12-17

References Repose®

Cushions and Mattress Overlay

- REP-6201100 Repose® Cushen (45x45 cm) + pump
- REP-6231100 Repose® Cushen (45x45 cm) + PU cover + pump
- REP-6101100 Repose® Mattresss Overlay (190x77 cm) + pump
- REP-6121100 Repose® Mattress Overlay + cushen + pump
- REP-6601100 Repose® Trolley Mattress Overlay (178x55 cm) + pump

Foot Protectors and Wedges

- REP-6501100 Repose® Foot protector (38x20 cm) (2 pcs. in 1 pump)
- REP-6021100 Repose® Foot protector large (51x28 cm) (2 pcs. in 1 pump)
- REP-6711100 Repose® Wedge (69x45 cm) + pump
- ORT-REP/HEEL/SPALK Heel splint for Repose® foot Protector

Pediatric Range

- REP-6401100 Repose® Babytherm (76x53 cm) + pomp
- REP-6801100 Repose® Pediatric Overlay Mattress (129x67 cm) + pump
- REP-6801101 Repose® Babynest (58x33x10 cm) + pump

Repose® - Total Care Pack

- REP-6882994 10 Repose® Mattress Overlays 3 Repose® Wedges
5 Repose® Cushens 1 big Repose® pump
4 Repose® Foot protectors 1 little Repose® pump

Repose® - Bulk products without pump

- REP-6111100 1 box of 10 Repose® Mattress Overlays + 1 pump
- REP-6200100 1 box of 10 Repose® Cushens + 1 pump
- REP-6500200 1 box of 20 Repose® Foot protectors (38 x 20 cm) + 1 pump
- REP-6500201 1 box of 10 Repose® Foot protectors L (51 x 28 cm) + 1 pump
- REP-6700102 1 box of 10 Repose® Wedges + 1 pump

Repose® Inbo Mattress

- REP-INBO-80 Repose® Inbo Mattress 197 x 80 x 14 cm
- REP-INBO-90 Repose® Inbo Mattress 197 x 90 x 14 cm



Repose®


An effective aid in the prevention and treatment of pressure ulcers and for pain relief

ENG - 2009/Repose

JMV
 produtos hospitalares

Centro Empresarial Sintra - Estoril VI, Edifício T1
 Estrada de Albarraque - Linhó
 2710-297 Sintra
 Tel: +351 21 924 82 90 Fax: +351 21 924 49 99
 info@jmv.com.pt
 www.jmv.com.pt

Visite-nos no nosso site www.jmv.com.pt



Repose® Cushion
 Inflated dimensions:
 L 450 mm x B 450 mm
 WL: 139 kg
 Ref.: REP-6201100

Repose® Mattress Overlay & Cushen
 Inflated dimensions:
 L 1.900 mm x B 770 mm
 L 450 mm x B 450 mm
 WL: 139 kg
 Ref.: REP-6121100

Repose® Mattress Overlay
 Inflated dimensions:
 L 1.900 mm x B 770 mm
 WL: 139 kg
 Ref.: REP-6101100

Repose® Trolley Mattress Overlay
 Inflated dimensions:
 L 1.900 mm x B 550 mm
 WL: 139 kg
 Ref.: REP-6601100

Repose® Babytherm
 Inflated dimensions:
 L 760 mm x B 530 mm
 WL: 3 kg
 Ref.: REP-6401100

Repose® Pediatric Mattress Overlay
 Inflated dimensions:
 L 1.290 mm x B 670 mm
 WL: 139 kg
 Ref.: REP-6801100

Repose® Foot Protector
 Inflated dimensions:
 • L 380 mm x B 200 mm x H 180 mm
 REP-6501100
 • L 510 mm x B 280 mm x H 180 mm
 REP-6021100
 WL: 139 kg

Repose® Babynest
 Inflated dimensions:
 L 680 mm x B 330 mm x H 100 mm
 WL: 3 kg
 Ref.: REP-6801101

Repose® Wedge
 Inflated dimensions:
 L 690 mm x B 450 mm
 WL: 139 kg
 Ref.: REP-6711100

WL = Recommended maximum weight limit

Repose® Total Care Pack

A completely new concept for your ward!
 It is essential to prevent pressure ulcers. To prevent means to apply the right product at the least suspicion or in case of non-bleachable redness.
 This is now possible with the Repose® Total Care Pack.
 The box contains:
 • 10 Repose® overlay mattresses
 • 5 Repose® cushions
 • 4 pair of Repose® foot protectors
 • 3 Repose® wedges
 • 1 large Repose® pump
 • 1 small Repose® pump




Repose® Inbo mattress

- The Repose® Inbo Mattress combines the advantages of the Repose® mattress overlay with those of a traditional decubitus prevention mattress.
- The integrated Repose® mattress overlay composed of 15 interconnected air cells ensure a permanently low interface pressure and minimized shear and friction forces.
- The rigid mattress frame allows patients to get safely out of the bed.
- The Repose® Inbo Mattress provides effective pressure redistribution for all persons at risk of developing pressure ulcers and for persons with pressure related tissue damage.
- An external valve allows a quick and safe inflation of the integrated Repose® mattress overlay.
- The included pump guarantees at any time the optimum internal pressure.
- A high-quality, bi-elastic incontinence cover ensures efficient cleaning and longevity.



Repose® : Your first aid against decubitus!
Static air filled in top-quality, 50µ thin polyurethane!
 Repose®, a total concept for your battle against decubitus. Repose® combines static air with an extremely elastic polyurethane film. One single air cell provides optimal pressure redistribution and minimization of shear forces across the total body contact area. Even small movements result in reciprocal beneficial interface pressure changes elsewhere. As a result, this thin, elastic PU-film forms in combination with air an autodynamic surface with excellent clinical evidence, proven by several clinical studies and numerous case studies. Repose® can be used for the prevention of pressure ulcers at all risks levels including very high risk patients and for the treatment of pressure related tissue damage with clinical supervision. Thanks to its easy handling and immediate availability, Repose® is the ideal alternative for a premature use of expensive dynamic devices.

- R** evolutionary total concept
- E** fficient in the prevention and treatment of pressure sores
- P** ermanent maximum pressure redistribution
- O** ne of a kind product
- S** afe and easy handling
- E** conomy-price



Our concept:
A good basic mattress + Repose®

Your advantage:

- Decrease in prevalence and incidence of pressure sores
- Optimal protection of the patient in bed and in the chair
- Optimal care for “all” patients
- Dramatic reduction in cost
- Time-saving for nursing staff
- Space-saving and low cost storage and transportation